



SAFETY DATA SHEET

page 1/5

VINTAGE

PRIME PRESS

Printing date: May 10, 2023

SECTION 1. Identification of the substance/mixture and of the company/undertaking

- 1.1 Product identifier
Trade Name:
VINTAGE "PRIME PRESS"
- 1.2 Relevant identified uses of the substance or mixture and uses advised against
Relevant identified uses: Dental material
Uses advised against: No further data
- 1.3 Details of the supplier of the safety data sheet
Company/Undertaking identification
Manufacturer's Name: SHOFU DENTAL GmbH
Address: An der Pönt 70, 40885 Ratingen, Germany
Phone: +49 (0) 2102-8664-0
Fax: +49 (0) 2102-8664-64
E-Mail: info@shofu.de
Section in Charge: Quality Management & Regulatory Affairs
- 1.4 Emergency Telephone Number
+49-2102-8664-53 (SHOFU DENTAL GmbH) 24 hours / 7 days

SECTION 2. Hazards identification

- 2.1 Classification of the substance or mixture
CLASSIFICATION (EC 1272/2008)
Not applicable.
- 2.2 Label elements
LABEL IN ACCORDANCE WITH (EC) NO.1272/2008
Void
HAZARD-DETERMINING COMPONENTS OF LABELLING
Void
SIGNAL WORD
Void
HAZARD STATEMENTS
Void
PRECAUTIONARY STATEMENTS
Void
- 2.3 Other hazards
Results of PBT and vPvB assessment
PBT: Not applicable.
vPvB: Not applicable.

SECTION 3. Composition/information on ingredients

- 3.1 Chemical characterization: Mixtures
- 3.2 Description: Mixture of substances listed below with nonhazardous additions.
- 3.3 Dangerous components: Void
Lithium silicate glass

(Contd. on page2)

Version Number 1

Date of issue: May 9, 2023



SAFETY DATA SHEET

page 2/5
VINTAGE
PRIME PRESS
Printing date: May 10, 2023

3.4 Additional information: For the wording of the listed risk phrases refer to section 2.

(Contd. of page1)

SECTION 4. First-aid measures

4.1 Description of first aid measures

(Grinding dust)

Eye contact: Rinse cautiously with water for several minutes. If eye irritation persists, get medical advice/attention.

Skin contact: Rinse with water.

Ingestion: Rinse mouth and seek medical advice if necessary.

Inhalation: Remove victim to fresh air and keep at rest in a position comfortable for breathing.

4.2 Most important symptoms and effects, both acute and delayed

No further relevant information available.

4.3 Indication of any immediate medical attention and special treatment needed

No further relevant information available.

SECTION 5. Fire-fighting measures

5.1 Extinguishing Media:

The product is not flammable.

5.2 Special hazards arising from the substance or mixture:

No further relevant information available.

5.3 Advice for firefighters:

No further relevant information available.

SECTION 6. Accidental release measures

6.1 Personal precautions, protective equipment and emergency procedures:

Not required.

6.2 Environmental Precautions:

No special measures required.

6.3 Methods and material for containment and cleaning Up:

Pick up mechanically.

6.4 Reference to other section:

See Section 7 for information on safe handling.

See Section 8 for information on personal protection equipment.

See Section 13 for disposal information.

SECTION 7. Handling and storage

7.1 Precautions for safe handling:

Use in accordance with instruction. In case of grinding and polishing , wear dust mask and safety glasses, and using dust respirator and local exhaust.

7.2 Conditions for safe storage, including any incompatibilities:

Store away from sunlight.

7.3 Specific end use(s):

No further relevant information available.

(Contd. on page3)

Version Number 1

Date of issue: May 9, 2023



SAFETY DATA SHEET

page 3/5
VINTAGE
PRIME PRESS
Printing date: May 10, 2023

(Contd. of page2)

SECTION 8. Exposure controls/personal protection

- 8.1 Control parameters:
Ingredients with limit values that require monitoring at the workplace:
The product does not contain any relevant quantities of materials with critical values that have to be monitored at the workplace.
- 8.2 Exposure controls:
Respiratory Protection:
Dust mask
- Skin Protection: Hand Protection
Not required.
- Eye Protection: Safety goggles

SECTION 9. Physical and chemical properties

- 9.1 Information on basic physical and chemical properties:
- | | |
|---|-----------------|
| Appearance/Odor/Colour: | Odorless ingot. |
| Odour threshold: | Not determined. |
| pH: | Not determined. |
| Melting point/freezing point: | Not determined. |
| Boiling point: | Not determined. |
| Flash point: | Not determined. |
| Evaporation rate: | Not determined. |
| Flammability (solid, gas): | Not applicable. |
| Upper/lower flammability or explosive limits: | Not determined. |
| Vapour pressure: | Not determined. |
| Vapour density: | Not determined. |
| Relative density: | Not determined. |
| Solubility: water solubility | Insoluble |
| Partition coefficient (n-octanol/water): | Not determined. |
| Auto-ignition temperature: | Not determined. |
| Decomposition temperature: | Not determined. |
| Viscosity: | Not determined. |
| Explosive properties: | Not applicable. |
| Oxidising properties: | Not applicable. |
- 9.2 Other information:
No further relevant information available.

SECTION 10. Stability and reactivity

- 10.1 Reactivity:
No further relevant information available.
- 10.2 Chemical stability:
Stable under normal temperatures and pressures.
- 10.3 Possibility of hazardous reactions:
No dangerous reactions known.

(Contd. on page4)



SAFETY DATA SHEET

page 4/5
VINTAGE
PRIME PRESS
Printing date: May 10, 2023

(Contd. of page3)

10.4 Condition to Avoid:

No further relevant information available.

10.5 Incompatible materials:

No further relevant information available.

10.6 Hazardous Decomposition Products:

None under normal conditions of storage and use.

SECTION 11. Toxicological information

11.1 Information on hazard classes as defined in Regulation (EC) No 1272/2008:

Acute toxicity: Based on available data, the classification criteria are not met.

Skin corrosion/irritation: Based on available data, the classification criteria are not met.

Eye damage/irritation: Based on available data, the classification criteria are not met.

Sensitization to the respiratory tract:

Based on available data, the classification criteria are not met.

Skin sensitization: Based on available data, the classification criteria are not met.

Germ cell mutagenicity/Genotoxicity:

Based on available data, the classification criteria are not met.

Carcinogenicity: Based on available data, the classification criteria are not met.

Reproductive toxicity: Based on available data, the classification criteria are not met.

Effects on or via lactation: Lack of data.

Specific target organ toxicity (single exposure):

Based on available data, the classification criteria are not met.

Specific target organ toxicity (repeated exposure):

Based on available data, the classification criteria are not met.

Aspiration hazard: Based on available data, the classification criteria are not met.

11.2 Information on other hazards:

No further relevant information available.

SECTION 12. Ecological information

12.1 Toxicity:

No further relevant information available.

12.2 Persistence and degradability:

No further relevant information available.

12.3 Bioaccumulative potential:

No further relevant information available.

12.4 Mobility in soil:

No further relevant information available.

12.5 Results of PBT and vPvB assessment:

Not applicable.

12.6 Endocrine disrupting properties:

This product does not contain any known or suspected endocrine disruptors.

12.7 Other adverse effects:

No further relevant information available.

(Contd. on page5)



SAFETY DATA SHEET

page 5/5
VINTAGE
PRIME PRESS
Printing date: May 10, 2023

(Contd. of page4)

SECTION 13. Disposal considerations

13.1 Waste treatment methods:

Dispose of contents/container to in accordance with local/regional/national/international regulations.

SECTION 14. Transport information

- | | |
|---|--|
| 14.1 UN number or ID number: | Void |
| 14.2 UN proper shipping name: | Void |
| 14.3 Transport hazard class (es): | Void |
| 14.4 Packing group: | Void |
| 14.5 Environmental hazards: | No further relevant information available. |
| 14.6 Special precautions for user: | Not applicable. |
| 14.7 Maritime transport in bulk according to IMO instruments: | Not applicable. |

SECTION 15. Regulatory information

- 15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture:
- EU REGULATIONS: See Section 2
 - Other regulations, limitations and prohibitive regulations:
The product is a medical device according to the EC-directive 93/42/EEC.
- 15.2 Chemical safety assessment:
A Chemical Safety Assessment has not been carried out.

SECTION 16. Other information

This product is intended for use by dental professionals (instrument/material).

Abbreviations and acronyms:

PBT: Persistent, Bioaccumulative and Toxic
vPvB: very Persistent and very Bioaccumulative