

Safety Data Sheet

according to Regulation (EC) No 1907/2006

ceraMotion Ti

Revision date: 15.04.2024

Product code: 644

Page 1 of 7

SECTION 1: Identification of the substance/mixture and of the company/undertaking**1.1. Product identifier**

ceraMotion Ti

Further trade names

Base Dentin, Base Dentin Mod, Bonder Powder, Correction, Dentin, Dentin Mod Chroma, Dentin Mod Fluo, Incisal, Incisal Mod, Incisal Opal, Incisal Transpa, Opaque, Opaque Mod, Shoulder, Transpa, Chroma Concept Opaque, Chroma Concept Dentin, Chroma Concept Incisal, Gingival, Glaze, Touch Up Base Dentin, Touch Up Dentin, Touch Up Incisal

Article No.: 251-021-20, 251-021-40, 251-022-20, 251-022-40, 251-023-20, 251-023-40, 251-024-20, 251-024-40, 251-025-20, 251-025-40, 251-026-20, 251-026-40, 251-027-20, 251-027-40, 251-028-20, 251-028-40, 251-029-20, 251-029-40, 251-030-20, 251-030-40, 251-031-20, 251-031-40, 251-032-20, 251-032-40, 251-033-20, 251-033-40, 251-034-20, 251-034-40, 251-035-20, 251-035-40, 251-036-20, 251-036-40, 251-107-20, 251-102-20, 251-104-20, 251-105-20, 251-103-20, 251-101-20, 251-106-20, 251-550-20, 251-550-40, 251-095-03, 251-041-20, 251-041-40, 251-042-20, 251-042-40, 251-043-20, 251-043-40, 251-044-20, 251-044-40, 251-045-20, 251-045-40, 251-046-20, 251-046-40, 251-047-20, 251-047-40, 251-048-20, 251-048-40, 251-049-20, 251-049-40, 251-050-20, 251-050-40, 251-051-20, 251-051-40, 251-052-20, 251-052-40, 251-053-20, 251-053-40, 251-054-20, 251-054-40, 251-055-20, 251-055-40, 251-056-20, 251-056-40, 251-111-20, 251-112-20, 251-113-20, 251-114-20, 251-115-20, 251-117-20, 251-116-20, 251-060-20, 251-060-40, 251-061-20, 251-061-40, 251-062-20, 251-062-40, 251-124-20, 251-123-20, 251-125-20, 251-121-20, 251-122-20, 251-065-20, 251-065-40, 251-066-20, 251-066-40, 251-067-20, 251-067-40, 251-070-20, 251-070-40, 251-071-20, 251-071-40, 251-072-20, 251-072-40, 251-001-20, 251-001-40, 251-002-20, 251-002-40, 251-003-20, 251-003-40, 251-004-20, 251-004-40, 251-005-20, 251-005-40, 251-006-20, 251-006-40, 251-007-20, 251-007-40, 251-008-20, 251-008-40, 251-009-20, 251-009-40, 251-010-20, 251-010-40, 251-011-20, 251-011-40, 251-012-20, 251-012-40, 251-013-20, 251-013-40, 251-014-20, 251-014-40, 251-015-20, 251-015-40, 251-016-20, 251-016-40, 251-081-20, 251-082-20, 251-083-20, 251-084-20, 251-086-20, 251-085-20, 251-075-20, 251-075-40, 251-018-20, 251-019-20, 251-221-20, 251-222-20, 251-223-20, 251-224-20, 251-231-20, 251-232-20, 251-233-20, 251-234-20, 251-235-20, 251-236-20, 251-237-20, 251-238-20, 251-239-20, 251-240-20, 251-241-20, 251-242-20, 251-245-20, 251-090-20, 251-091-20, 251-092-20, 251-093-20, 251-251-20, 251-252-20, 251-253-20, 251-255-20, 251-256-20, 251-257-20, 251-260-20, 251-261-20, 251-262-20, 251-270-02

1.2. Relevant identified uses of the substance or mixture and uses advised against**Use of the substance/mixture**

Health services.

1.3. Details of the supplier of the safety data sheet

Company name:	DENTAURUM GmbH & Co. KG		
Street:	Turnstr. 31		
Place:	D-75228 Ispringen		
Telephone:	+49 7231 803 0	Telefax: +49 7231 803 295	
e-mail:	info@dentaurum.de		
Contact person:	Chemie		
Internet:	www.dentaurum.com		

1.4. Emergency telephone number: +497231803-184 7:00-16:15 (Mo-Do) 7:00-13:15 (Fr)
+497231803-0 16:15-18:00 (Mo-Do) 13:15-18:00 (Fr)

SECTION 2: Hazards identification**2.1. Classification of the substance or mixture**

No risks worthy of mention.

2.2. Label elements

Safety Data Sheet

according to Regulation (EC) No 1907/2006

ceraMotion Ti

Revision date: 15.04.2024

Product code: 644

Page 2 of 7

2.3. Other hazards

No information available.

SECTION 3: Composition/information on ingredients**3.2. Mixtures****Chemical characterization**

Synthetic vitroceraamic powder with following ingredients:

SiO₂, Al₂O₃, K₂O, Na₂O, Li₂O, CaO, BaO, MgO, B₂O₃, CaF₂, Ca₂P₂O₇, SnO₂, ZrO₂, BaTiO₃ + Pigments**Hazardous components**

none (according to Regulation (EC) No 1907/2006 (REACH))

SECTION 4: First aid measures**4.1. Description of first aid measures****After inhalation**

Provide fresh air. Get medical advice/attention if you feel unwell.

After contact with skin

After contact with skin, wash immediately with: Water and soap.

After contact with eyes

In case of contact with eyes flush immediately with plenty of flowing water for 10 to 15 minutes holding eyelids apart and consult an ophthalmologist.

After ingestion

Rinse mouth immediately and drink 1 glass of water.

4.2. Most important symptoms and effects, both acute and delayed

No information available.

4.3. Indication of any immediate medical attention and special treatment needed

Treat symptomatically.

SECTION 5: Firefighting measures**5.1. Extinguishing media****Suitable extinguishing media**

Co-ordinate fire-fighting measures to the fire surroundings.

5.2. Special hazards arising from the substance or mixture

The product itself does not burn.

5.3. Advice for firefighters

Use suitable breathing apparatus.

SECTION 6: Accidental release measures**6.1. Personal precautions, protective equipment and emergency procedures****General advice**

Avoid dust formation.

6.2. Environmental precautions

No special environmental measures are necessary. Clean contaminated articles and floor according to the environmental legislation.

6.3. Methods and material for containment and cleaning up

ceraMotion Ti

Revision date: 15.04.2024

Product code: 644

Page 3 of 7

Other information

Take up mechanically.

6.4. Reference to other sections

Safe handling: see section 7

Personal protection equipment: see section 8

Disposal: see section 13

SECTION 7: Handling and storage**7.1. Precautions for safe handling****Advice on safe handling**

Generation/formation of dust,

Advice on general occupational hygiene

Do not breathe dust. Wash hands before breaks and after work.

7.2. Conditions for safe storage, including any incompatibilities**Requirements for storage rooms and vessels**

Keep/Store only in original container.

Further information on storage conditions

Keep container tightly closed. Store in a cool dry place.

7.3. Specific end use(s)

Health services.

SECTION 8: Exposure controls/personal protection**8.1. Control parameters****8.2. Exposure controls****Individual protection measures, such as personal protective equipment****Eye/face protection**

Wear eye/face protection.

Hand protection

Wear protective gloves.

Recommendation:

Dermatril Single-use gloves. NBR (Nitrile rubber). 0,11 mm

When handling with chemical substances, protective gloves must be worn with the CE-label including the four control digits. The quality of the protective gloves resistant to chemicals must be chosen as a function of the specific working place concentration and quantity of hazardous substances. For special purposes, it is recommended to check the resistance to chemicals of the protective gloves mentioned above together with the supplier of these gloves.

Skin protection

Use of protective clothing.

Respiratory protection

Use appropriate respiratory protection.

SECTION 9: Physical and chemical properties**9.1. Information on basic physical and chemical properties**

Physical state:

Powder

Colour:

various, depending on coloration

Odour:

odourless

Melting point/freezing point:

approx. 860 °C

ceraMotion Ti

Revision date: 15.04.2024

Product code: 644

Page 4 of 7

Boiling point or initial boiling point and boiling range:	not determined
Flammability:	not applicable
Lower explosion limits:	not determined
Upper explosion limits:	not determined
Flash point:	not determined
Auto-ignition temperature:	not determined
Decomposition temperature:	not determined
pH-Value:	not applicable
Water solubility:	The study does not need to be conducted because the substance is known to be insoluble in water.
Solubility in other solvents	
not determined	
Partition coefficient n-octanol/water:	not determined
Vapour pressure:	not determined
Density:	approx. 2,4 - 2,7 g/cm ² g/cm ³
Relative vapour density:	not determined

9.2. Other information

Information with regard to physical hazard classes

Explosive properties

The product is not: Explosive.

Oxidizing properties

The product is not: oxidising.

Other safety characteristics

Evaporation rate:

not determined

Solid content:

not determined

SECTION 10: Stability and reactivity

10.1. Reactivity

No hazardous reaction when handled and stored according to provisions.

10.2. Chemical stability

The product is stable under storage at normal ambient temperatures.

10.3. Possibility of hazardous reactions

No risks worthy of mention.

10.4. Conditions to avoid

none

10.5. Incompatible materials

No information available.

10.6. Hazardous decomposition products

No known hazardous decomposition products.

Further information

Thermal decomposition: none

SECTION 11: Toxicological information

11.1. Information on hazard classes as defined in Regulation (EC) No 1272/2008

Acute toxicity

Based on available data, the classification criteria are not met.

Safety Data Sheet

according to Regulation (EC) No 1907/2006

ceraMotion Ti

Revision date: 15.04.2024

Product code: 644

Page 5 of 7

Irritation and corrosivity

Based on available data, the classification criteria are not met.

Sensitising effects

Based on available data, the classification criteria are not met.

Carcinogenic/mutagenic/toxic effects for reproduction

Based on available data, the classification criteria are not met.

STOT-single exposure

Based on available data, the classification criteria are not met.

STOT-repeated exposure

Based on available data, the classification criteria are not met.

Aspiration hazard

Based on available data, the classification criteria are not met.

11.2. Information on other hazards**Endocrine disrupting properties**

This product does not contain a substance that has endocrine disrupting properties with respect to humans as no components meets the criteria.

Further information

Toxicological data are not available.

SECTION 12: Ecological information**12.1. Toxicity**

The product is not: Ecotoxic.

12.2. Persistence and degradability

The product has not been tested.

12.3. Bioaccumulative potential

The product has not been tested.

12.4. Mobility in soil

The product has not been tested.

12.5. Results of PBT and vPvB assessment

The substances in the mixture do not meet the PBT/vPvB criteria according to REACH, annex XIII.

The product has not been tested.

12.6. Endocrine disrupting properties

This product does not contain a substance that has endocrine disrupting properties with respect to non-target organisms as no components meets the criteria.

12.7. Other adverse effects

Avoid release to the environment.

Further information

No information available.

SECTION 13: Disposal considerations**13.1. Waste treatment methods****Disposal recommendations**

Dispose of waste according to applicable legislation. May be disposed of in household waste landfill.

List of Wastes Code - contaminated packaging

200301 MUNICIPAL WASTES (HOUSEHOLD WASTE AND SIMILAR COMMERCIAL, INDUSTRIAL AND INSTITUTIONAL WASTES) INCLUDING SEPARATELY COLLECTED FRACTIONS; other municipal wastes; mixed municipal waste

Safety Data Sheet

according to Regulation (EC) No 1907/2006

ceraMotion Ti

Revision date: 15.04.2024

Product code: 644

Page 6 of 7

Contaminated packaging

Non-contaminated packages may be recycled.

SECTION 14: Transport information**Land transport (ADR/RID)**

<u>14.1. UN number or ID number:</u>	No dangerous good in sense of this transport regulation.
<u>14.2. UN proper shipping name:</u>	No dangerous good in sense of this transport regulation.
<u>14.3. Transport hazard class(es):</u>	No dangerous good in sense of this transport regulation.
<u>14.4. Packing group:</u>	No dangerous good in sense of this transport regulation.

Inland waterways transport (ADN)

<u>14.1. UN number or ID number:</u>	No dangerous good in sense of this transport regulation.
<u>14.2. UN proper shipping name:</u>	No dangerous good in sense of this transport regulation.
<u>14.3. Transport hazard class(es):</u>	No dangerous good in sense of this transport regulation.
<u>14.4. Packing group:</u>	No dangerous good in sense of this transport regulation.

Marine transport (IMDG)

<u>14.1. UN number or ID number:</u>	No dangerous good in sense of this transport regulation.
<u>14.2. UN proper shipping name:</u>	No dangerous good in sense of this transport regulation.
<u>14.3. Transport hazard class(es):</u>	No dangerous good in sense of this transport regulation.
<u>14.4. Packing group:</u>	No dangerous good in sense of this transport regulation.

Air transport (ICAO-TI/IATA-DGR)

<u>14.1. UN number or ID number:</u>	No dangerous good in sense of this transport regulation.
<u>14.2. UN proper shipping name:</u>	No dangerous good in sense of this transport regulation.
<u>14.3. Transport hazard class(es):</u>	No dangerous good in sense of this transport regulation.
<u>14.4. Packing group:</u>	No dangerous good in sense of this transport regulation.

14.5. Environmental hazards

ENVIRONMENTALLY HAZARDOUS: No

14.6. Special precautions for user

No information available.

14.7. Maritime transport in bulk according to IMO instruments

not applicable

Other applicable information

No dangerous good in sense of these transport regulations.

SECTION 15: Regulatory information**15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture****EU regulatory information****Additional information**

For this product, a hazard label is not required according to section 1.3.4 of Annex I of the CLP regulation.

National regulatory information

Water hazard class (D): - - non-hazardous to water

15.2. Chemical safety assessment

Chemical safety assessments for substances in this mixture were not carried out.

SECTION 16: Other information**Changes**

This data sheet contains changes from the previous version in section(s):
1,2,3,4,5,6,7,8,9,10,11,12,13,14,15,16.

Safety Data Sheet

according to Regulation (EC) No 1907/2006

ceraMotion Ti

Revision date: 15.04.2024

Product code: 644

Page 7 of 7

Abbreviations and acronyms

CLP: Classification, labelling and Packaging
REACH: Registration, Evaluation and Authorization of Chemicals
GHS: Globally Harmonised System of Classification, Labelling and Packaging of Chemicals
UN: United Nations
CAS: Chemical Abstracts Service
DNEL: Derived No Effect Level
DMEL: Derived Minimal Effect Level
PNEC: Predicted No Effect Concentration
ATE: Acute toxicity estimate
LC50: Lethal concentration, 50%
LD50: Lethal dose, 50%
LL50: Lethal loading, 50%
EL50: Effect loading, 50%
EC50: Effective Concentration 50%
ErC50: Effective Concentration 50%, growth rate
NOEC: No Observed Effect Concentration
BCF: Bio-concentration factor
PBT: persistent, bioaccumulative, toxic
vPvB: very persistent, very bioaccumulative
ADR: Accord européen sur le transport des marchandises dangereuses par Route
(European Agreement concerning the International Carriage of Dangerous Goods by Road)
RID: Regulations concerning the international carriage of dangerous goods by rail
ADN: European Agreement concerning the International Carriage of Dangerous Goods by Inland Waterways
(Accord européen relatif au transport international des marchandises dangereuses par voies de navigation intérieures)
IMDG: International Maritime Code for Dangerous Goods
EmS: Emergency Schedules
MFAG: Medical First Aid Guide
IATA: International Air Transport Association
ICAO: International Civil Aviation Organization
MARPOL: International Convention for the Prevention of Marine Pollution from Ships
IBC: Intermediate Bulk Container
SVHC: Substance of Very High Concern
For abbreviations and acronyms, see: ECHA Guidance on information requirements and chemical safety assessment, chapter R.20 (Table of terms and abbreviations).

Further Information

The above information describes exclusively the safety requirements of the product and is based on our present-day knowledge. The information is intended to give you advice about the safe handling of the product named in this safety data sheet, for storage, processing, transport and disposal. The information cannot be transferred to other products. In the case of mixing the product with other products or in the case of processing, the information on this safety data sheet is not necessarily valid for the new made-up material.

(The data for the hazardous ingredients were taken respectively from the last version of the sub-contractor's safety data sheet.)